





Article

Development and Validation of Competitive ELISA for Detection of H5 Hemagglutinin Antibodies

Orie Hochman¹, Wanhong Xu¹ , Ming Yang¹, Chengbo Yang² , Aruna Ambagala¹ , Anna Rogiewicz², Joseph J. Wang³ and Yohannes Berhane^{1,2,4,*} 

¹ National Centre for Foreign Animal Disease, 1015 Arlington Street, Winnipeg, MB R3E 3M4, Canada; orie.hochman@inspection.gc.ca (O.H.); wanhong.xu@inspection.gc.ca (W.X.); ming.yang@inspection.gc.ca (M.Y.); aruna.ambagala@inspection.gc.ca (A.A.)

² Department of Animal Science, University of Manitoba, Winnipeg, MB R3T 2N2, Canada; chengbo.yang@umanitoba.ca (C.Y.); anna.rogiewicz@umanitoba.ca (A.R.)

³ BioStone Animal Health, Southlake, TX 76092, USA; jwang@biostoneah.com

⁴ Department of Veterinary Pathology, Western College of Veterinary Medicine, University of Saskatchewan, Saskatoon, SK S7N 5B4, Canada

* Correspondence: yohannes.berhane@inspection.gc.ca; Tel.: +1-204-789-7062

Abstract: Influenza A viruses (IAVs) belonging to the goose/Guangdong (Gs/GD)-lineage H5Nx remain a major concern for the global poultry industry, wildlife, and humans. The hemagglutinin (HA) protein is the dominant antigenic epitope carrier within IAV, which in turn triggers substantial immunogenic responses in the infected host. The current study describes the development and validation of a highly sensitive competitive H5 ELISA (cELISA) based on a novel monoclonal antibody developed in mice immunized with inactivated virus H5N1 (A/Turkey/ON/6213/1966). The cELISA is capable of detecting the H5 antibody response to a wider range of H5-subtype viruses belonging to both North American and Eurasia lineages, including the Gs/GD H5Nx from clade 2.3.4.4b that is currently causing the highly pathogenic avian influenza outbreaks in Eurasia, Africa, and Latin and North America. The developed H5 cELISA provides a specific, sensitive, and species-independent serological assay for the rapid detection of H5 antibodies. The assay is more robust and more sensitive than the hemagglutination inhibition assay, which is the “Gold standard”. The assay can be used in serological diagnosis, serosurveillance, and vaccine monitoring of serum samples collected from different species of animals.

Keywords: avian influenza; competitive ELISA; hemagglutinin; H5 subtype; serodiagnosis



Citation: Hochman, O.; Xu, W.; Yang, M.; Yang, C.; Ambagala, A.; Rogiewicz, A.; Wang, J.J.; Berhane, Y. Development and Validation of Competitive ELISA for Detection of H5 Hemagglutinin Antibodies. *Poultry* **2023**, *2*, 349–362. <https://doi.org/10.3390/poultry2030026>

Academic Editors: Priscilla F. Gerber and Surya Paudel

Received: 6 May 2023

Revised: 29 June 2023

Accepted: 30 June 2023

Published: 7 July 2023



Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<https://creativecommons.org/licenses/by/4.0/>).

1. Introduction

Influenza A virus is classified as being part of the family of *Orthomyxoviridae* and the genus of Influenza A virus, and it contains a negative-sense, single-stranded RNA genome composed of eight gene segments encoding at least 11 viral proteins [1,2]. IAVs are divided into low-pathogenic (LP) and highly pathogenic (HP) pathotypes based on the genetic features and/or disease severity in chickens. To date, 16 hemagglutinin (HA) and 9 neuraminidase (NA) subtypes of IAVs have been isolated from aquatic and domestic birds [3]. However, only the H5 and H7 subtypes of IAV have been found to possess the potential to mutate into highly pathogenic viruses during circulation in poultry [4–6]. All H5 and H7 IA viruses have been classified as notifiable viruses by the World Organization for Animal Health (WOAH) since 2005 [7]. Highly pathogenic avian influenza viruses (HP-IAVs) that originated from low-pathogenic avian influenza viruses (LPIAVs) have mostly been eradicated from poultry via stamping-out programs, except for the Gs/GD lineage of H5, which rose to prominence in 1996. The Gs/GD lineage of H5 viruses has spread and been maintained in poultry and wild aquatic bird reservoirs until today, producing infections in poultry, wild birds, and mammals, including humans in many countries in

Asia, Africa, Europe, and North America [8,9]. Currently, H5 HPIAVs have been classified into clades 0–9 and different subclades based on phylogenetic analysis of their HA genes [10]. The clade 2.3.4.4 HA of H5 viruses has evolved into eight subclades, namely, clades 2.3.4.4a to 2.3.4.4h [11]. Since the H5N1 viruses bearing the clade 2.3.4.4b HA gene emerged in October 2020 in the Netherlands, the viruses have spread to many countries in Europe, Africa, Asia, and North America and are responsible for the loss of millions of domestic poultry [12–15]. The same viruses emerged in Canada in the winter of 2021 through migratory wild birds and have spread to all 10 provinces, resulting in a loss of over 7 million domestic birds [16,17].

To prevent the extensive circulation of notifiable IAVs in both wild birds and domestic poultry populations, mandated serological surveillance and survey programs have been implemented worldwide. Classical laboratory serologic tools include agar gel immunodiffusion (AGID), the hemagglutination inhibition (HI) test, and the virus neutralization (VN) test. The AGID test is used to detect group-specific antigens of the influenza A virus (i.e., nucleoprotein and matrix proteins). It is a simple and economical assay but it is time-consuming and not suitable for large-scale specimen screening [18]. The HI assay is considered to be the reference test by WOA. It is the subtype-specific test for influenza serologic differential diagnosis, but the downsides of this test include antigen selection and titration, nonspecific inhibitor removal, and a long processing time [19]. The VN test is recommended for the identification of HPIAV using live infectious viruses. Thus, the test needs to be performed in higher-level containment labs and has significant limitations in rapid and high-throughput diagnostics [19]. Conversely, enzyme-linked immunosorbent assay (ELISA) is widely accepted as an essential method that is superior in regard to throughput, turnaround time, and accuracy for early diagnosis [20]. Among the 11 proteins encoded by the IAV, haemagglutinin (HA) and nucleoprotein (NP) are particularly targeted for serological diagnosis due to their antigenic abundance [21]. An indirect ELISA is considered a valuable and inexpensive test routinely used for the detection of antibodies to IAV in birds. However, this technique is not feasible when screening samples from a wide range of bird species. Alternatively, the competitive ELISA (cELISA) is suitable for broader use in multiple animal species as a species-specific enzyme conjugate is not required [22].

The success of a reliable cELISA requires a highly specific monoclonal antibody (mAb) to recognize a broadly conserved and dominant antigenic epitope. In addition to the availability of two commercial multispecies NP-based cELISA kits (i.e., IDVet, Montpellier, France, and Idexx, Westbrook, ME, USA), a recombinant NP-based cELISA developed in our laboratory is routinely used as a screening tool for IAV serological monitoring owing to the group-specific feature of the NP antigen [23]. Previously, Yang and colleagues in our laboratory successfully generated a high-affinity, subtype-specific mAb #10 (F37H5N1-56) based on immunizing mice with the inactivated American-H5 (clade Am-nonGs/GD) LPIAV strain, A/Turkey/ON/6213/1966 (H5N1). This mAb has been shown to react with sera raised from a wide range of HP and LP H5 strains of North American and Eurasian origins [24]. In this study, a cELISA immunoassay based on mAb #10 (F37H5N1-56) and the recombinant H5 HA antigen derived from a Eurasian strain A/Teal/Germany/Wv632/2005 (H5N1) (clade EA-nonGs/GD) was developed, validated and evaluated with high specificity, sensitivity, and low variability. In particular, the developed cELISA is capable of detecting antibodies in serum samples collected from different species of animals that were infected with currently circulating clade 2.3.4.4b H5Nx viruses.

2. Materials and Methods

2.1. Preparation of Viruses

IAVs of the H5 subtype used in this study were acquired from the National Centre for Foreign Animal Disease (NCFAD) viral depository. Allantoic fluid from virus-inoculated embryonated chicken eggs was collected and tested by hemagglutination assay (HA) as previously described [24]. All procedures were performed under biocontainment level (BSL) 3 conditions.

2.2. Preparation and Expression of Recombinant HA Proteins

The preparation and expression of recombinant HA protein were described previously [25]. For the development of the competitive ELISA, we used a recombinantly expressed, full-length hemagglutinin protein of A/Teal/Germany/Wv632/2005 (H5N1) (GenBank accession no. CY061885). This full-length H5 hemagglutinin gene was cloned into a pAB- beeTM-FH vector (AB Vector, LLC, San Diego, CA, USA) by GenScript (GenScript USA Inc., Piscataway, NJ, USA). The vector containing the desired gene sequences was then purified and co-transfected with linearized baculovirus vector DNA, ProFoldTM-ER1 (AB Vector, LLC, San Diego, CA, USA), into *Spodoptera frugiperda* (Sf9) insect cells to generate a recombinant baculovirus containing a full-length HA gene. This recombinant baculovirus was plaque-purified and sequenced to verify the correct nucleotide sequence. The titer of the recombinant H5 baculoviral stock (PFU/mL) was determined by a viral plaque assay. *Trichoplusia ni* (Tni) insect cells were then infected with the recombinant virus at an MOI of 10. The infected insect cells were incubated for 72 h at 27 °C via shaking and were harvested via centrifugation (10,000 rpm for 20 min). The cell pellet was lysed using the I-PER insect cell protein extraction reagent (Pierce Biotechnology, Rockford, IL, USA) on ice for 10 min, and then was centrifuged again at 15,000 rpm for 20 min. The soluble protein in the supernatant was purified using Ni-NTA resin. The recombinant protein was detected and confirmed via Western blotting using a 6× histidine-specific mAb (Novagen, Billerica, MA, USA). After the determination of the correct protein size, the concentration of the antigen was determined using the Qubit 2.0 Flex and bicinchoninic acid assay (BCA). An average between the two methods was used to calculate the antigen-per-well concentration.

2.3. Hemagglutination Inhibition (HI) Assay

The HI assay was performed according to the OIE manual [26]. The HI titer was expressed as the highest dilution of serum, resulting in the complete inhibition of 4 hemagglutinin units (HAUs) of the virus. The HI titer was regarded as being positive if there was inhibition at a serum dilution of 1:16 or more against 4 HAUs of the virus. Hemagglutination inhibition titers for the H5-positive chicken antisera were obtained using a homologous virus.

2.4. Sera Used in This Study

Chickens and turkeys of specific-pathogen-free flocks were used in the production of Influenza A virus antisera. All procedures involving experimental animal inoculations and care complied with the Canadian Council of Animal Care guidelines. These animals were allowed to acclimate for 5–7 days prior to inoculations. For LPIA viruses, fresh allantoic fluid was collected that contained an HA of greater than 1:16, which was used and diluted 1:10 in PBS. For HPAI viruses, fresh allantoic fluid with an HA of greater than 1:16 was collected and inactivated via binary ethylenimine inactivation (BEI) or β-propiolactone (BPL) virus inactivation procedures, as stated above. The inactivated HPAI viruses were then formulated with 50 µg of QuilA/mL of allantoic fluid prior to inoculation. Before inoculation of the birds, the undiluted virus was plated on blood agar plates and incubated at 37 °C ± 2 °C overnight. If plates presented positive for bacterial presence, the inoculum was then filtered through a 0.45 µm syringe filter. If not already diluted, the inoculum was diluted 1:10 in PBS prior to inoculation (LPIA viruses were already previously diluted).

LPIA viruses were inoculated intravenously via the wing vein in conjunction with three drops of inoculum that were administered via the oral, nasal, and ocular routes of the birds. HPAI viruses were injected intramuscularly along with Emulsigen adjuvant into different regions of the pectoral muscle. Approximately 0.25 mL of HPAI inoculum was administered into each injection site. A predetermined schedule was made in consultation with veterinarians for bleeding. These serum samples were then collected from each infected/immunized bird and tested for the presence/absence of antibodies against the inoculated virus by HI assay. Booster doses were determined using the HI reactivity results. Once the desired antibody response and titer had been achieved, the animals

were anesthetized and exsanguinated via predesignated protocols. Serum pools were then created and tested again for HI reactivity to obtain a greater than 1:16 positive result. The serum was then aliquoted and stored at -20°C for later usage.

Sera used to assess the diagnostic sensitivity and specificity include experimental anti-H5 sera and field sera submitted through the Canadian Notifiable Avian Influenza Surveillance System (CanNAISS). The IAV-positive status was confirmed by IAV-NP cELISA [23] and the subtype status was confirmed by HI assay.

Sera used for the kinetic time course antibody response study were derived from 4 specific-pathogen-free (SPF) leghorn chickens that were inoculated with 0.5 mL of the β -propiolactone-inactivated antigen (A/Teal/Germany/Wv632/2005 (H5N1)) by intramuscular injection and housed in the same room. Chicken sera were collected at 0 days post-vaccination (DPI), 14 DPI, 21 DPI, 28 DPI, 35 DPI, and 42 DPI.

Reference chicken sera that represent 16 serotypes of avian influenza virus (IAV); 3 paramyxoviruses (APMV), types 1, 2, and 3; and avian metapneumovirus subtypes A, B, and C were used to assess cross-reactivity in the rec-H5 cELISA. The virus strains used to produce 16 IAV reference sera were A/Swine/Iowa/31 (H1N1), A/Pintail/AB/293/77 (H2N9), A/Duck/UK/1/63 (H3N8), A/Duck/Czech/56 (H4N6), A/Turkey/BC/FAV5-52/09 (H5N2), A/Turkey/Mass/65 (H6N2), A/Chicken/BC/514/04 (H7N3), A/Turkey/ON/6118/162 (H8N4), A/Turkey/MN/12877/1285/81 (H9N2), A/Chicken/Germany/N/49 (H10N7), A/Duck/Memphis/546/74 (H11N9), A/Duck/AB/60/76 (H12N5), A/Gull/MD/704/77/BEL-42 (H13N1), A/Mallard/Gurjev/263/82 (H14N5), A/Duck/Australia/341/83 (H15N8) and A/BHG/Sweden/5/99 (H16N3). The virus strains used to produce three APMV reference sera were Newcastle disease virus (NDV), LaSota (AMPV-1), Yucaipa (AMPV-2), and Ty6661/67 (AMPV-3). The virus strains used to produce three aMPV reference sera were 14/1 UK (aMPV-A), Hungary/657/4 (aMPV-B), and Colorado (aMPV-C).

2.5. Selection of Anti-H5 HA Monoclonal Antibody

A panel of nine H5-subtype monoclonal antibodies generated from mice immunized with the inactivated virus A/Turkey/ON/6213/66 (H5N1) was described previously [24]. The characteristics of the mAbs are summarized as follows: Eight out of nine had conformational epitopes, including mAb #10 (F37H5N1-56). The #10 mAb was one of the two mAbs that reacted with all thirteen H5 strains tested in the indirect ELISA, double-antibody sandwich (DAS) ELISA, and dot blot assay. Eight out of nine mAbs were non-neutralizing mAbs by HI assay, including mAb #10. mAb #10 was one of the three mAbs that worked in immunohistochemistry and immunofluorescence assays for H5 antigen detection. The isotype of mAb #10 is IgG2b/ κ . Given the characteristics of mAb #10 described above, it was chosen to explore its suitability within the cELISA immunoassay format in this study.

2.6. Competitive ELISA Condition Parameters

The cELISA condition parameters were described previously [25]. The optimum conditions of the cELISA were assessed by various factors, such as microtiter plate type, blocking/diluent buffer, pre-assay blocking, incubation temperatures and/or agitation, coating antigen/mAb, and testing serum concentration.

The determination of the optimal amount of coating antigen per well was achieved by utilizing a checkerboard-style titration by indirect ELISA. The amount of coating antigen per well is reported in $\mu\text{g/mL}$, and was calculated from the stock protein concentration that was determined by finding the average between the Qubit 2.0 and BCA methods. The concentration of monoclonal and secondary antibodies was determined by a checkerboard titration and expressed as $\mu\text{g/mL}$ per well.

A two-fold dilution series of neat positive H5 IAV serum and negative normal chicken serum at a ratio ranging from a 1:2.5 to 1:40 was undertaken and compared.

Two different variants of a 96-well microtiter plate were used to determine the best platform for the cELISA test. A non-treated, flat-bottom, clear polystyrene plate (Thermo Fisher Scientific, Rochester, NY, USA) that has a naturally hydrophobic surface property

was tested in parallel with a specially treated, flat-bottom polystyrene plate that was designed to have a high-protein-binding surface property for maximum capture (Maxisorp, Thermo Fisher). Assays with the plate types were run in parallel using 14 reference sera in duplicate to determine differences in the resulting outcomes using the identical protocol and lot-numbered reagents.

Utilizing a blocking buffer as a pre-step treatment was studied to determine its effects on the assay. This procedure was carried out on non-treated 96-well microtiter plates with 3% heat-inactivated FBS in PBS-T as a blocking buffer.

Blocking/diluent buffers tested included 3% heat-inactivated fetal bovine serum (FBS) in PBS-T, 5% skim milk in PBS-T, and 2% bovine + 2% rabbit serum in PBS-T. These buffers were tested on both antigens against the same seven H5 chicken antisera in duplicate to observe any changes the buffers may cause in the assay.

The effect of agitation during the incubation periods of the cELISA was examined both at 37 °C and at an ambient temperature.

The effect of the antigen coating time was assessed to determine if the overnight coating was necessary due to its time restraint in regard to testing delays and planning, or if coating the antigen at a higher temperature for a shorter time would yield similar results. In this regard, antigens coated overnight at 4 °C and antigens coated for 2 h at 37 °C were tested and compared.

2.7. Development of H5 cELISA

The working concentration of the recombinant H5 proteins was diluted in a 0.06 M carbonate buffer with pH 9.6 (45 mM sodium bicarbonate and 18 mM sodium carbonate). Each well of non-treated 96-well microtiter plates was coated overnight at 4 °C with 100 µL of the above-diluted antigen and then was washed 5× with PBS plus 0.05% Tween20 (PBS-T). Each well was blocked with 100 µL of 3% heat-inactivated FBS in PBS-T and incubated for 1 h at 37 °C. After washing 5× with PBS-T, equal volumes (50 µL) of test serum samples and competitive mAb diluted to the working concentration with the blocking buffer were mixed and incubated simultaneously at 37 °C for 1 h. Wells containing competitive mAb mixed with standard anti-H5 sera, negative sera, or no serum were utilized as a positive, negative, or dilution control (DC), respectively. Each serum sample was tested in duplicate. After 5× washing, 100 µL of horseradish peroxidase-labeled goat anti-mouse IgG conjugate (Jackson ImmunoResearch Laboratories, West Grove, PA, USA) was added at a dilution of 1/2000, incubated for 1 h at 37 °C, and then washed and developed with 3,3',5,5'-tetramethylbenzidine as a peroxidase substrate (TMB) (Sigma-Aldrich, St. Louis, MO, USA). The reaction was carried out for 15 min at room temperature, which was followed by adding 50 µL of 2 M H₂SO₄ to each well to stop the reaction. Colorimetric development was quantified spectrophotometrically at 450 nm with a Molecular Devices EMax precision microplate reader (Molecular Devices, LLC, San Jose, CA, USA). Results were interpreted as the percentage of inhibition (PI) calculated according to the following formula: % inhibition = $((1 - (\text{OD}_{450} \text{ of serum sample} / \text{OD}_{450} \text{ of DC})) \times 100$.

2.8. Validation Methods for the Determination of Diagnostic Sensitivity, Specificity, Repeatability, and Threshold Cut-Off Level

To assess the diagnostic sensitivity and specificity of the assay, all sera with known serostatus (confirmed by HI assay) were used. This included experimental chicken and turkey serum samples produced at the NCFAD and field serum samples submitted to NCFAD through CanNAISS. For validating the assay, confirmed negative serum samples and positive serum samples were used in a receiver operating characteristic (ROC) analysis and calculated to assess diagnostic performance, which included the determination of sensitivity, specificity, and the threshold cut-off using SigmaPlot version 14.0 (Systat Software, San Jose, CA, USA).

The repeatability of the assay was assessed by running the reference sera in multiple replicates within the same run or between runs. The intra-assay repeatability was calcu-

lated for 45 replicates on 2 separate plates and then repeated over 2 days for inter-assay repeatability assessment. The values were expressed as a mean, standard deviation, and percent coefficient of variation (CV%) for the repeated measure.

2.9. Statistical Analysis

The rec-H5 antigen-based cELISA was validated using ROC analysis and was used to calculate diagnostic performance, which included the determination of sensitivity, specificity, and the threshold cut-off using SigmaPlot version 14.0 (Systat Software, San Jose, CA, USA) for statistical analysis. The measurement of agreement between the gold-standard HI test and the cELISA tests was assessed using Cohen's kappa (κ) test. To interpret the given κ value, the criteria laid out by Landis and Koch was utilized [27]. $\kappa \leq 0.00$ was designated as poor agreement, $0.00 < \kappa \leq 0.20$ as slight agreement, $0.21 < \kappa \leq 0.40$ as fair agreement, $0.41 < \kappa \leq 0.60$ as moderate agreement, $0.61 < \kappa \leq 0.80$ as substantial agreement, and $0.81 < \kappa < 1.00$ as almost perfect agreement.

3. Results

3.1. Expression of Recombinant Influenza A Virus Subtype H5 Hemagglutinin Protein

The rec-H5 protein expressed in *Tni* insect cells was harvested and purified, and the expression of the H5 protein was confirmed with Western blotting using a 6× histidine-specific mAb (Figure 1A) as well as an anti-H5 reference chicken serum (Figure 1B). Due to the glycosylation of the proteins during expression in insect cells, the molecular weight of the protein was slightly different from the calculated weight (69 kDa).




Figure 1. Western blot analysis showing the detection of recombinantly expressed H5 HA protein: (A) detected with anti-polyhistidine mAb; (B) detected with reference serum produced from A/Chicken/PA/1370/1/83 (H5N2). MW, molecular weight; 1, recombinant H5 HA protein derived from HA gene of A/Teal/Germany/Wv632/2005 (H5N1) generated in this study; 2, reference recombinant H5 HA protein derived from HA gene of A/Canada goose/OR/AH0012452/2015 (H5N8); 3, uninfected cell lysates.

3.2. H5 cELISA Optimization

To optimize the H5 cELISA, various parameters were assessed. The amount of coating antigen per well was assessed using a titration curve (Figure 2A). The dilution series did not significantly improve or impair the percent of inhibition of the assay but significantly affected the optical density. The appropriate concentration of antigen used for coating

was determined at 58 ng/well due to the higher percent inhibition (PI) and optimal OD value when this concentration was used. The concentrations of mAb #10 and the secondary antibody were determined by checkerboard titrations. The optimal concentration for the mAb was 28 ng/well, whereas, for the anti-mouse-HRP conjugate, it was 43 ng/well. To determine the proper working dilution of serum for efficient assay development, the serial dilution of a positive H5 IAV serum and negative normal chicken serum was conducted. As shown in Figure 2B, the antisera diluted at 1:10 in dilution buffer provided the optimal condition, with the lowest negative serum background interference and highest PI value. In addition, we compared differences in the type of microtiter plates used, types of blocking/diluent buffers, blocking vs. non-blocking conditions, shaking vs. stationary incubation, and antigen coating overnight at 4 °C versus 2 h at 37 °C on their impacts on PI value (data not shown). Optimal cELISA results were obtained using non-treated, flat-bottom, and clear polystyrene plates coated overnight at 4 °C with 3% heat-inactivated FBS in PBST as the blocking/diluent buffer and with constant agitation during the incubation.




Figure 2. Antigen and antiserum optimization for the H5 cELISA. (A) Coating antigen curve comparing antigen concentration to optical density (OD) and percent inhibition (PI) values. Monoclonal antibody F37 #10 and anti-H5 chicken serum from chickens infected with A/Chicken/PA/1370/1/83 were used at a fixed concentration. The final coating antigen concentration was based on an OD value ranging between 1.0 and 2.0 along with the highest PI value. (B) Optimal antisera working dilution for testing was determined via 2-fold serial dilution from 1:2.5 to 1:40. The working dilution condition that provided the highest PI values on positive H5 IAV serum along with the lowest PI values on normal chicken serum was determined to be 1:10 dilution. The black dotted vertical line indicates the working dilution condition. The rec-H5 antigen was used at 58 ng/well.

3.3. The Diagnostic Sensitivity, Specificity, and Threshold Cut-Off Level

To determine the sensitivity, specificity, and threshold cut-off level, a ROC analysis was performed using antisera from multiple avian species (sera from naturally or experimentally infected chickens, turkeys, and wild birds) of known serostatus. The influenza A antibody status was first evaluated by IAV-NP ELISA, and all positive sera were further tested with HI assay to determine if the sera were H5 positive. A total of 275 experimentally produced IAV-H5-positive sera and 1204 IAV-NP-negative sera were used for the validation of the H5 cELISA with a rec-H5 full-length protein derived from the A/Teal/Germany/Wv632/2005 (H5N1) virus as the coating antigen. The optimal cut-off value and corresponding sensitivity and specificity of the rec-H5 IAV cELISA are presented in Figure 3. The rec-H5 cELISA demonstrated a sensitivity of 99.6% and a specificity of 97.6% at a cut-off value of 35%. Hence, a sample run on the rec-H5 cELISA with a PI% value of greater or equal to 35% was considered positive.




Figure 3. Receiver operator characteristic (ROC) validation and determination of diagnostic sensitivity and specificity of the rec-H5 cELISA. (A) Dot histogram pair for the rec-H5 cELISA using 275 IAV-H5-positive sera and 1204 IAV-NP-negative sera. The horizontal line bisecting the dot plots represents the cut-off value that gives the optimal diagnostic sensitivity and specificity. (B) Graph of sensitivity and specificity vs. cut-off using data generated from ROC analysis. The arrow indicating the cut-off was set at the intersection of the sensitivity and specificity plots, where optimal diagnostic sensitivity and specificity were achieved at 99.6% and 97.6%, respectively.

3.4. Analytic Sensitivity of rec-H5 cELISA

The analytic sensitivity of the rec-H5 cELISA was evaluated using three reference anti-H5 chicken sera in comparison to the gold-standard HI assay (Figure 4). The detection limit of the rec-H5 cELISA was below the diagnostic detection limit of the HI assay (2^4 HI titer) when tested in comparison with three known H5-positive chicken sera. The detection limit of the rec-H5 cELISA was the 2^0 HI titer, while the 2^4 HI titer was considered a positive result in the HI test. Hemagglutination inhibition titers for the H5-positive chicken antisera used to determine analytical sensitivity were obtained using a homologous virus.




Figure 4. Analytic sensitivity of the rec-H5 cELISA compared to HI assay. Positive chicken sera raised from A/Teal/Germany/Wv632/05, A/Chicken/PA/1370/83, or A/Turkey/MN/3689-1551/81 were two-fold diluted and subjected to cELISA and HI assay for H5 antibody detection. The horizontal dotted lines indicate the cut-off value of rec-H5 cELISA.

3.5. Assessment of Repeatability

The repeatability of our rec-H5 cELISA was assessed by testing reference sera in duplicate within the same run or between runs. The intra-assay repeatability was calculated using 45 replicates on two separate plates, and then was repeated over 2 days for inter-assay repeatability assessment. The rec-H5 cELISA showed a percent coefficient of variation (CV%) with 7.1% intra-assay variability and 5.5% inter-assay variability. The results demonstrate that the tests are highly repeatable for diagnostic applications.

3.6. Reactivity of the mAb F37 #10 and rec-H5 Antigen to Different H5 Strains and Assessment for Cross-Reactivity to Other Avian-Based Viruses

The broadness of the developed rec-H5 cELISA was assessed against diverse H5 strains (Table 1). The rec-H5 cELISA can detect H5-subtype antibodies collected from different species of birds experimentally infected with both North American- and Eurasian-lineage H5-subtype viruses including the currently circulating Gs/GD lineage H5Nx viruses belonging to the 2.3.4.4b clade.

Table 1. Summary of rec-H5 cELISA detection.

Reference Serum Strain Name	Subtype	Anti-Serum Source	Lineage/Clade	Pathogenicity	* PI (%)
A/Mallard/Nova Scotia/FAV-102-3/2022	H5N1	Duck	Gs/GD 2.3.4.4b	HPAI	71
A/Red-breasted Merganser/British Columbia/FAV-899-13/2022	H5N1	Duck	Gs/GD 2.3.4.4b	HPAI	73
A/chicken/Newfoundland and Labrador/FAV-33/2021	H5N1	Goose	Gs/GD 2.3.4.4b	HPAI	90
A/Chicken/Vietnam/14/2005	H5N1	Chicken	Gs/GD 2.3.2	HPAI	70
A/Turkey/British Columbia/FAV2/2015	H5N1	Goose	Gs/GD 2.3.4.4c	HPAI	89
A/Turkey/Ontario/6213/1966	H5N1	Chicken	Am_nonGs/GD	LPAI	85
A/Swan/Germany/R065/2006	H5N1	Chicken	Gs/GD 2.2	HPAI	97
A/Turkey/British Columbia/FAV10/2014	H5N2	Chicken	Gs/GD 2.3.4.4c	HPAI	96
A/Duck/British Columbia/CN26-6/2005	H5N2	Chicken	Am_nonGs/GD	LPAI	96
A/Quail/Oregon/20719/1986	H5N2	Chicken	Am_nonGs/GD	LPAI	83
A/Chicken/Pennsylvania/1370/1983	H5N2	Chicken	Am_nonGs/GD	HPAI	94
A/Turkey/Minnesota/3689-1551/1981	H5N2	Chicken	Am_nonGs/GD	LPAI	66
A/Chicken/Washington/13413/1984	H5N2	Chicken	Am_nonGs/GD	HPAI	84
A/Turkey/California/21860/1984	H5N3	Chicken	Am_nonGs/GD	LPAI	80
A/Tern/SouthAfrica/1961	H5N3	Chicken	EA_nonGs/GD	HPAI	87
A/Turkey/Ireland/1378/1983	H5N8	Chicken	EA_nonGs/GD	HPAI	89
A/Turkey/Wisconsin/1/1968	H5N9	Chicken	Am_nonGs/GD	LPAI	73

*, Percentage of inhibition.

The potential for the rec-H5 cELISA to cross-react with antibodies developed against other avian virus sera was examined. Representative antiserum with subtypes ranging from IAV H1 to H16, avian metapneumovirus (aMPV) subtypes A, B and C, and paramyxovirus (APMV) types 1, 2 and 3 were used to determine cross-reactivity. As shown in Figure 5, none of the IAV-, APMV- and aMPV-positive sera cross-reacted in the rec-H5 cELISA based on a cut-off value of 35% except for the control H5 antiserum.




Figure 5. Assessment of potential cross-reactivity with antibodies developed against various non-H5 IAV, one LPIAV H5, or other avian viruses. Antiserum from representative viruses was described in the Material and Methods section. The horizontal dotted line indicates a cut-off PI (30%) for IAV NP cELISA and the horizontal solid line represents a cut-off PI (35%) for rec-H5 cELISA.

3.7. The Kinetics of Anti-H5 Antibody Response

In Figure 6, we show the kinetics of the mean anti-H5 antibodies detected by the H5 cELISA during the follow-up after the infection period (42 days). The early positive H5-specific antibody response was detected at 7 dpi. The anti-H5 antibodies peaked at 21 dpi and remained positive until the end of the experiment on day 42.




Figure 6. Kinetic time course antibody evaluation. Antibody time course kinetics were expressed as the mean percentage of inhibition in the H5 cELISA using serum samples from chickens experimentally infected with A/Teal/Germany/Wv632/2005 (H5N1). The dashed horizontal line indicates the diagnostic cut-off for the test. The chickens were boosted at 28 DPI.

3.8. Diagnostic Performance of the rec-H5 cELISA Using Field Samples

To assess the capability of our rec-H5 cELISA to detect H5 antibodies in wild birds naturally infected with H5-subtype viruses, a total of 48 serum samples collected from wild birds in 2022 were first screened by an in-house-developed IAV NP cELISA. Forty-seven out of forty-eight samples tested positive for NP cELISA. All forty-eight samples were then tested via H5 cELISA and HI assay. For the HI test, the 2.3.4.4b virus, A/Turkey/ON/FAV-0162-144/2022 and A/Teal/Germany/Wv632/2005 (H5N1) were used as antigens. Consistent with the NP cELISA results, the same 47 samples tested positive using A/Turkey/ON/FAV-0162-144/2022 (H5N1) as an antigen, but none of the samples tested positive in the HI assay using A/Teal/Germany/Wv632/2005 (H5N1). This indicates that the ducks were potentially exposed to Gs/Gd lineage H5-subtype viruses from the 2.3.4.4b clade. The agreement between the H5 cELISA and HI assay was assessed by Cohen's κ analysis, and a perfect agreement was observed between the two methods ($\kappa = 1$).

4. Discussion

Influenza A viruses of the H5 subtype are widely distributed and have caused numerous outbreaks in poultry and wild birds worldwide [28]. Highly pathogenic avian influenza viruses (HPIAVs) of the H5 subtype are a constant threat to public health due to their potential to cause zoonotic disease. Avian influenza outbreaks associated with Gs/Gd-lineage H5Nx viruses belonging to clade 2.3.4.4b have been reported in Europe, Africa, Asia and America since 2020 and have continued to threaten the poultry industry and wildlife, accounting for an unprecedentedly large amount of mortality in birds [17]. Effective prevention and control of H5 outbreaks require active serological surveillance in wild birds and poultry. Since 2008, Canada has designed and implemented the Canadian Notifiable Avian Influenza Surveillance System (CanNAISS) with six surveillance activities (i.e., wild bird surveillance; passive surveillance in domestic poultry when clinical signs suggestive of notifiable avian influenza are reported; targeted surveillance when notifiable avian influenza is detected; pre-slaughter surveillance in commercial poultry; hatchery supply flock surveillance; and voluntary enhanced surveillance in the poultry genetic exporters sector) aiming to detect, prevent and eliminate the presence of H5 and H7 subtypes (notifiable avian influenza) in Canada's domestic poultry flocks [28]. Due to the increase in the sampling for serosurveillance, we would like to replace the labor-intensive HI assay with a simple ELISA assay for determining serotype specificity that is more cost-effective and more robust, and that does not require the use of live IAVs as in the HI assay.

In the current study, we describe the development and validation of the H5 cELISA which was developed using a recombinant Eurasian-origin, full-length H5 protein as a coating antigen and a competitor H5 mAb developed from mice immunized with the North American-lineage H5N1 strain. Using the H5 cELISA, we were able to detect H5 antibodies developed against different H5 strains belonging to both North American- and Eurasian-lineage viruses including Gs/GD-lineage H5-subtype viruses belonging to 2.3.4.4c and 2.3.4.4b clades. We have, by far, extensively expanded our detection spectrum compared to previously described assays [24,29–31]. The H5 mAb used in this study was previously characterized as an isotype of the IgG2b/ κ that binds with a conformational epitope of the H5 HA protein [24]. In this study, we further confirmed that mAb #10 does not demonstrate anti-hemagglutinin activity and is not suitable for use in the HI assay. We confirmed that the epitope location of this mAb is within the HA1 subunit of the HA protein, as it did not react with the recombinant HA2 protein in a dot blot or indirect ELISA test (data not shown). The absence of anti-hemagglutinin activity or features in this mAb might have contributed to its broader reactivity spectrum, which is possibly due to a smaller burden of immune selection pressure on its binding epitope. Studies have shown that frequent residue alterations at the HA antigenic sites were observed with vaccinations implemented in poultry [32,33]. Thus, if the mAb epitope is within the antigenic sites, it has more chance to lose reactivity due to antigenic evolution.

The Gaussian distribution approach is commonly used to determine the cut-off value for diagnostic tests, in which a cut-off value is defined as the mean plus two standard deviations of the negative reference samples. This approach ensures a diagnostic specificity without considering the diagnostic sensitivity [34]. In this study, we determined the cut-off value by use of ROC analysis to allow for a combined measure of diagnostic sensitivity and specificity. The HI is regarded as the gold-standard subtyping test for IAV in routine diagnostics recommended by the WOAHA. The H5 cELISA developed in this study showed higher accordance with the HI assay for H5 antibody detection with the determined cut-off value.

The newly developed cELISA has proven to be more sensitive compared to the HI assay and also more suitable for the mass screening of serum samples collected from different species of animals. In the hemagglutinin inhibition assay, a serum sample with an HI titer ≥ 8 is considered positive. The detection limit of the H5 cELISA was determined using three serially diluted H5-positive chicken sera that were tested with the HI assay against homologous H5-subtype viruses. The detection limit of the assay was determined to be below the diagnostic detection limit of the HI assay for all three serum samples tested, indicating that the cELISA was superior in analytical sensitivity compared to the HI assay. Despite being the “gold standard” assay for subtyping IAVs, the HI assay has some drawbacks compared to the ELISA. These include: (1) the HI assay is more laborious and time-consuming in selecting the virus antigen, red blood cell preparations, serum dilution and the manual reading of results; (2) the selection of a homologous viral antigen is difficult when the infection history of the animals from which the serum samples were collected is unknown; and (3) due to the zoonotic potential of IAVs, all procedures involving live virus antigens require the use of a biological safety cabinet. In contrast to the HI assay, the developed H5 cELISA in our study is easy to automate for the large-scale monitoring of H5-subtype-specific antibodies in sera samples collected from different species of animals. In addition, the established rec-H5 cELISA has a biosafety advantage over the HI assay as it can be performed in a regular lab without the use of a biological safety cabinet.

An attractive feature of the developed H5 cELISA is its ability to detect H5 antibodies of the currently circulating H5Nx viruses belonging to clade 2.3.4.4b. To assess the agreement between the developed H5 cELISA and HI assay, we tested 48 field serum samples collected from mallard ducks in 2022 using both assays. A perfect agreement was observed in Cohen's κ analysis. Experimental studies have shown that mallard ducks can excrete abundant viruses without clinical or pathologic evidence of debilitating disease, supporting the potential role of these ducks as long-distance vectors of HPIAV [35]. When these migratory birds are in close proximity to domestic poultry facilities, poultry birds can become infected through direct contact with secretions or feces of infected birds, contaminated surfaces, or infected food and water supplies. The H5 subtype responsible for the outbreak that began in 2021 is highly pathogenic and thus is of heightened concern. The ability of our H5 cELISA to detect anti-H5 antibodies in the sera derived from the current H5N1 HPIAVs makes it a rapid and effective serological screening tool for the early detection and control of subclinical H5 infection. The H5 cELISA could be used in the serosurveillance of H5 antibodies in different species of birds and also in mammalian species. The assay could also be used for monitoring the vaccination status of domestic poultry where the disease is endemic and vaccination is used to control the spread of the virus.

5. Conclusions

This study demonstrated that our rec-H5 cELISA can be used as an alternative to the HI assay for serologic diagnosis of serum samples collected from animals infected with H5-subtype viruses. This rec-H5 cELISA has an inherently wider capability to detect H5 antibodies in serum samples collected from different species of animals pre-exposed to both North American- and Eurasian-lineage H5-subtype viruses. The rapid turnaround time and higher throughput ability coupled with its highest specificity, sensitivity, and

reproducibility mean this assay is a promising approach for massive epidemiological screening and the serosurveillance of antibody response against the H5 subtype viruses.

Author Contributions: Conceptualization, W.X., Y.B. and M.Y.; methodology, W.X.; validation, O.H.; formal analysis, O.H. and W.X.; investigation, O.H. and W.X.; resources, Y.B. and M.Y.; data curation, O.H. and W.X.; writing—original draft preparation, O.H. and W.X.; writing—review and editing, W.X., M.Y., C.Y., A.A., A.R., J.J.W. and Y.B.; visualization, O.H., W.X., M.Y., C.Y., A.A., A.R., J.J.W. and Y.B.; supervision, Y.B., M.Y. and C.Y.; project administration, Y.B.; funding acquisition, Y.B. All authors have read and agreed to the published version of the manuscript.

Funding: This research was funded by the Technology Development of the Canadian Food Inspection Agency, project number WIN-A-1814, and O.H. was funded by BioStone Animal Health for pursuing a Master of Science degree (BioStone Animal Health LLC (#01970)).

Institutional Review Board Statement: The original samples and animal experiments were processed at the National Centre for Foreign Animal Disease laboratories based on the Canadian Science Centre for Human & Animal Health standard operating procedures under ISO 17025. All animal procedures were approved by the Canadian Science Centre for Human and Animal Health Animal Care committee under AUD C-18-001 on 12 January 2018.

Informed Consent Statement: Not applicable.

Data Availability Statement: Not applicable.

Acknowledgments: The authors thank the technical personnel of NCFAD for providing animal care and husbandry during animal experiments. The authors thank Tamiru Alkie, Kathleen Hooper-McGrevy, and Charles Nfon for their review of the manuscript.

Conflicts of Interest: The authors declare no conflict of interest. The funders had no role in the design of the study; in the collection, analyses or interpretation of data; in the writing of the manuscript; or in the decision to publish the results.

References

1. Tollis, M.; Di Trani, L. Recent Developments in Avian Influenza Research: Epidemiology and Immunoprophylaxis. *Vet. J.* **2002**, *164*, 202–215. [[CrossRef](#)] [[PubMed](#)]
2. Cheung, T.K.W.; Poon, L.L.M. Biology of Influenza A Virus. *Ann. N. Y. Acad. Sci.* **2007**, *1102*, 1–25. [[CrossRef](#)] [[PubMed](#)]
3. Fouchier, R.A.M.; Munster, V.; Wallensten, A.; Bestebroer, T.M.; Herfst, S.; Smith, D.; Rimmelzwaan, G.F.; Olsen, B.; Osterhaus, A.D.M.E. Characterization of a Novel Influenza A Virus Hemagglutinin Subtype (H16) obtained from Black-Headed Gulls. *J. Virol.* **2005**, *79*, 2814–2822. [[CrossRef](#)] [[PubMed](#)]
4. Deshpande, K.L.; Fried, V.A.; Ando, M.; Webster, R.G. Glycosylation Affects Cleavage of an H5N2 Influenza Virus Hemagglutinin and Regulates Virulence. *Proc. Natl. Acad. Sci. USA* **1987**, *84*, 36–40. [[CrossRef](#)]
5. Capua, I.; Alexander, D.J. Avian Influenza Infections in Birds—A Moving Target. *Influenza Other Respir. Viruses* **2007**, *1*, 11–18. [[CrossRef](#)]
6. Bonfanti, L.; Monne, I.; Tamba, M.; Santucci, U.; Massi, P.; Patregnani, T.; Loli Piccolomini, L.; Natalini, S.; Ferri, G.; Cattoli, G.; et al. Highly Pathogenic H7N7 Avian Influenza in Italy. *Vet. Rec.* **2014**, *174*, 382. [[CrossRef](#)] [[PubMed](#)]
7. World Organization for Animal Health (OIE). Chapter 2.7.12 Avian Influenza. In *Terrestrial Animal Health Code*, 6th ed.; OIE: Paris, France, 2007; pp. 279–285.
8. Rohm, C.; Horimoto, T.; Kawaoka, Y.; Suss, J.; Webster, R.G. Do Hemagglutinin Genes of Highly Pathogenic Avian Influenza Viruses Constitute Unique Phylogenetic Lineages? *Virology* **1995**, *209*, 664–670. [[CrossRef](#)]
9. Swayne, D.E.; Suarez, D.L.; Sims, L. Influenza. In *Diseases of Poultry*; Swayne, D.E., Boulianne, M., Logue, C., McDougald, L.D., Nair, V., Suarez, D.L., Eds.; Wiley: Ames, IA, USA, 2020; pp. 210–256.
10. WHO/OIE/FAO H5N1 Evolution Working Group. Toward a Unified Nomenclature System for Highly Pathogenic Avian Influenza Virus (H5N1). *Emerg. Infect. Dis.* **2008**, *14*, e1. [[CrossRef](#)]
11. Li, Y.; Li, M.; Li, Y.; Tian, J.; Bai, X.; Yang, C.; Shi, J.; Ai, R.; Chen, W.; Zhang, W.; et al. Outbreaks of Highly Pathogenic Avian Influenza (H5N6) Virus Subclade 2.3.4.4h in Swans, Xinjiang, Western China, 2020. *Emerg. Infect. Dis.* **2020**, *26*, 2956–2960. [[CrossRef](#)]
12. Bevins, S.N.; Shriner, S.A.; Cumbee, J.C.J.; Dilione, K.E.; Douglass, K.E.; Ellis, J.W.; Killian, M.L.; Torchetti, M.K.; Lenocho, J.B. Intercontinental Movement of Highly Pathogenic Avian Influenza A(H5N1) Clade 2.3.4.4 Virus to the United States, 2021. *Emerg. Infect. Dis.* **2022**, *28*, 1006–1011. [[CrossRef](#)]

13. European Food Safety Authority; European Centre for Disease Prevention and Control; European Union Reference Laboratory for Avian Influenza; Adlhoch, C.; Fusaro, A.; Gonzales, J.L.; Kuiken, T.; Marangon, S.; Niqueux, E.; Staubach, C.; et al. Avian Influenza Overview December 2021–March 2022. *EFSA J.* **2022**, *20*, e07289. [[PubMed](#)]
14. Sagong, M.; Lee, Y.; Song, S.; Cha, R.M.; Lee, E.; Kang, Y.; Cho, H.; Kang, H.; Lee, Y.; Lee, K. Emergence of Clade 2.3.4.4b Novel Reassortant H5N1 High Pathogenicity Avian Influenza Virus in South Korea during Late 2021. *Transbound Emerg. Dis.* **2022**, *69*, e3255–e3260. [[CrossRef](#)] [[PubMed](#)]
15. Wille, M.; Barr, I.G. Resurgence of Avian Influenza Virus. *Science* **2022**, *376*, 459–460. [[CrossRef](#)]
16. Alkie, T.N.; Lopes, S.; Hisanaga, T.; Xu, W.; Suderman, M.; Koziuk, J.; Fisher, M.; Redford, T.; Lung, O.; Joseph, T.; et al. A Threat from both Sides: Multiple Introductions of Genetically Distinct H5 HPAI Viruses into Canada Via both East Asia–Australasia/Pacific and Atlantic Flyways. *Virus Evol.* **2022**, *8*, veac077. [[CrossRef](#)] [[PubMed](#)]
17. WOA. Available online: <https://wahis.Woah.Org/> (accessed on 23 April 2023).
18. Charlton, B.; Crossley, B.; Hietala, S. Conventional and Future Diagnostics for Avian Influenza. *Comp. Immunol. Microbiol. Infect. Dis.* **2009**, *32*, 341–350. [[CrossRef](#)]
19. Katz, J.M.; Hancock, K.; Xu, X. Serologic Assays for Influenza Surveillance, Diagnosis and Vaccine Evaluation. *Expert Rev. Anti Infect. Ther.* **2011**, *9*, 669–683. [[CrossRef](#)]
20. Okamatsu, M.; Hiono, T.; Kida, H.; Sakoda, Y. Recent Developments in the Diagnosis of Avian Influenza. *Vet. J.* **2016**, *215*, 82–86. [[CrossRef](#)]
21. Webster, R.G.; Bean, W.J.; Gorman, O.T.; Chambers, T.M.; Kawaoka, Y. Evolution and Ecology of IAV. *Microbiol. Rev.* **1992**, *56*, 152–179. [[CrossRef](#)]
22. Spackman, E.; Killian, M.L. Detection of Influenza a Antibodies in Avian Samples by ELISA. *Methods Mol. Biol.* **2020**, *2123*, 177–193.
23. Yang, M.; Berhane, Y.; Salo, T.; Li, M.; Hole, K.; Clavijo, A. Development and Application of Monoclonal Antibodies Against Avian Influenza Virus Nucleoprotein. *J. Virol. Methods* **2008**, *147*, 265–274. [[CrossRef](#)]
24. Yang, M.; Clavijo, A.; Graham, J.; Salo, T.; Hole, K.; Berhane, Y. Production and Diagnostic Application of Monoclonal Antibodies Against Influenza Virus H5. *J. Virol. Methods* **2009**, *162*, 194–202. [[CrossRef](#)] [[PubMed](#)]
25. Hochman, O. Development and Validation of a Recombinant-H5 Hemagglutinin-Based Competitive ELISA for Serodiagnosis of Avian Influenza A Subtype H5 Antibodies. Available online: <https://mspace.Lib.Umanitoba.ca/handle/1993/36905> (accessed on 1 April 2023).
26. World Organization for Animal Health (OIE). OIE Manual of Diagnostic Tests and Vaccines for Terrestrial Animals. In *Avian Influenza*; OIE: Paris, France, 2008; Volume 1, pp. 465–481.
27. Landis, J.R.; Koch, G.G. An Application of Hierarchical Kappa-Type Statistics in the Assessment of Majority Agreement among Multiple Observers. *Biometrics* **1977**, *33*, 363–374. [[CrossRef](#)] [[PubMed](#)]
28. Harfoot, R.; Webby, R.J. H5 Influenza, a Global Update. *J. Microbiol.* **2017**, *55*, 196–203. [[CrossRef](#)]
29. Dlugolenski, D.; Hauck, R.; Hogan, R.J.; Michel, F.; Mundt, E. Production of H5-Specific Monoclonal Antibodies and the Development of a Competitive Enzyme-Linked Immunosorbent Assay for Detection of H5 Antibodies in Multiple Species. *Avian Dis.* **2010**, *54*, 644–649. [[CrossRef](#)]
30. Postel, A.; Ziller, M.; Rudolf, M.; Letzel, T.; Ehrlich, R.; Pourquier, P.; Dauber, M.; Grund, C.; Beer, M.; Harder, T.C. Broad Spectrum Reactivity Versus Subtype Specificity-Trade-Offs in Serodiagnosis of Influenza A Virus Infections by Competitive ELISA. *J. Virol. Methods* **2011**, *173*, 49–59. [[CrossRef](#)] [[PubMed](#)]
31. Shehata, A.A.; Fiebig, P.; Sultan, H.; Hafez, M.; Liebert, U.G. Development of a Recombinant ELISA using Yeast (*Pichia Pastoris*)-Expressed Polypeptides for Detection of Antibodies Against Avian Influenza A Subtype H5. *J. Virol. Methods* **2012**, *180*, 18–25. [[CrossRef](#)] [[PubMed](#)]
32. Xu, W.; Navarro-Lopez, R.; Solis-Hernandez, M.; Liljehult-Fuentes, F.; Molina-Montiel, M.; Lagunas-Ayala, M.; Rocha-Martinez, M.; Ferrara-Tijera, E.; Perez de la Rosa, J.; Berhane, Y. Evolutionary Dynamics of Mexican Lineage H5N2 Avian Influenza Viruses. *Viruses* **2022**, *14*, 958. [[CrossRef](#)]
33. Navarro-Lopez, R.; Xu, W.; Gomez-Romero, N.; Velazquez-Salinas, L.; Berhane, Y. Phylogenetic Inference of the 2022 Highly Pathogenic H7N3 Avian Influenza Outbreak in Northern Mexico. *Pathogens* **2022**, *11*, 1284. [[CrossRef](#)]
34. Greiner, M.; Pfeiffer, D.; Smith, R.D. Principles and Practical Application of the Receiver-Operating Characteristic Analysis for Diagnostic Tests. *Prev. Vet. Med.* **2000**, *45*, 23–41. [[CrossRef](#)]
35. Keawcharoen, J.; van Riel, D.; van Amerongen, G.; Bestebroer, T.; Beyer, W.E.; van Lavieren, R.; Osterhaus, A.D.M.E.; Fouchier, R.A.M.; Kuiken, T. Wild Ducks as Long-Distance Vectors of Highly Pathogenic Avian Influenza Virus (H5N1). *Emerg. Infect. Dis.* **2008**, *14*, 600–607. [[CrossRef](#)]

Disclaimer/Publisher’s Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.