

# Supplementary Materials: A Toolbox of Genetically Encoded FRET-Based Biosensors for Rapid L-Lysine Analysis

Victoria Steffen, Julia Otten, Susann Engelmann, Andreas Radek, Michael Limberg, Bernd W. Koenig, Stephan Noack, Wolfgang Wiechert and Martina Pohl

## 1. Titration Curves of the Sensor Prototype with L-lysine, L-arginine, L-histidine, and L-glutamine




**Figure S1.** Binding isotherms of the sensor prototype with L-lysine, L-arginine, L-histidine, and L-glutamine. Distortion of the titration curve  $>0.1$  mM L-arginine and L-lysine, respectively, is due to the use of a crude sensor preparation. The results demonstrate that the sensor with the cpLAO-binding protein still has similar affinities for L-lysine ( $107 \mu\text{M}$ ) and L-arginine (ca.  $25 \mu\text{M}$ ), whereas the affinity for L-histidine is much lower ( $>100$  mM).

## 2. pH-Dependent Binding Isotherms of the Sensor Prototype




**Figure S2.** Binding isotherms of the sensor prototype in MOPS buffer with pH-values between pH 6.7 to 7.2. The solution of the sensor prototype was stored in 20 mM MOPS buffer, pH 7.3. This buffer was replaced by 20 mM MOPS buffer with the respective target pH via ultrafiltration. The recorded binding isotherms are shown in (a) and for better comparison the normalized binding isotherms are shown in (b).

### 3. Isothermal Calorimetry




**Figure S3.** Isothermal titration calorimetry data reflecting L-lysine binding to the binding protein cpLAO-BP (a) and the complete sensor prototype containing this binding protein (b). The upper two panels show baseline subtracted raw data for titration of 22  $\mu\text{M}$  cpLAO-BP with 100  $\mu\text{M}$  L-lysine (left) and 60  $\mu\text{M}$  sensor prototype with 1 mM L-lysine (right). The lower panels display enthalpy changes per added mole of L-lysine as a function of the total molar ratio of lysine to protein in the calorimeter cell, reflecting binding isotherms. Only the sigmoidal isotherm in the left panel can be reliably fit to a binding model (single binding site,  $K_d = 1.5 \mu\text{M}$ ,  $N = 0.36$ ,  $\Delta H = -58 \text{ kcal mol}^{-1}$ ). Please note the different scaling of the y-axis in (a) and (b).

### 4. Protein Sequences from *S. typhimurium* and *E. coli*

*S. typhimurium* ALPQTVRIGT DTTYAPFSSK DAKGEFI GFD IDLGNECMCKR MQVKCTWVAS DFDALIPSLK  
*E. coli* ALPETVRIGT DTTYAPFSSK DAKGDFVGF IDLGNECMCKR MQVKCTWVAS DFDALIPSLK

*S. typhimurium* AKKIDAISS LSITDKRQQE IAFSDKLYAA DSRLIAAKGS PIQPTLES LK GKHVGV LQGS  
*E. coli* AKKIDAISS LSITDKRQQE IAFSDKLYAA DSRLIAAKGS PIQPTLD SLK GKHVGV LQGS

*S. typhimurium* TQEAYANDNW RTKGVDVVAY ANQDLI YSDL TAGRLDAALQ DEVAASEGFL KQPAGKEYAF  
*E. coli* TQEAYANETW RSKGVDVVAY ANQDLVYSDL AAGRDLAALQ DEVAASEGFL KQPAGKDFAF

*S. typhimurium* AGPSVKDKKY FGDGTGVGLR KDDTELKAAF DKALTEL RQD GTYDKMAKKY FDFNVYGD  
*E. coli* AGSSVKDKKY FGDGTGVGLR KDDAELTA AF NKALGEL RQD GTYDKMAKKY FDFNVYGD

**Figure S4.** Comparison of the protein sequences derived from *S. typhimurium* and *E. coli*. The lysine binding proteins show high similarity. Homologous substitutions are marked in green, and non-homologous substitutions are marked in red. The amino acids involved in lysine binding are underlined. The sequence deleted during the circular permutation is marked in yellow.

### 5. Characteristic Mutations of the Fluorescent Proteins Relative to GFP

ECFP: F64L/S65T/Y66W/N146I/M153T/V163A [1]  
 Citrine: S65G/V68L/Q69M/S72A/T203Y/H231L [2]

## 6. DNA and Protein Sequence of the Sensor Construct with cpLAO-BP

The DNA sequence encoding the His-Tag is highlighted in **bold**, the sequence of CFP is marked blue, the recognition site of the restriction enzymes is underlined, the LAO-binding protein sequence is shown in green, and the Citrine sequence is shown in yellow. The same code was used for the protein sequence.

### 6.1. Nucleotide Sequence:

ATGCGGGGTTCT**CATCATCATCATCATCAT**GGTATGGCTGATACTCGCATTGGTGTAAC  
AATCTATAAGTCGGCTGGT**ATGGTGAGCAAGGGCGAGGAGCTGTTCACCGGGGTGGTGCC**  
**ATCCTGGTCGAGCTGGACGGCGACGTAAACGGCCACAAGTTCAGCGTGTCCGGCGAGGGC**  
**GAGGGCGATGCCACCTACGGCAAGCTGACCCTGAAGTTCATCTGCACCACCGGCAAGCTG**  
**CCCGTGCCCTGGCCCCACCCTCGTGACCACCCTGACCTGGGGCGTGCAGTGCTTCAGCCGCT**  
**ACCCCGACCACATGAAGCAGCACGACTTCTTCAAGTCCGCCATGCCCGAAGGCTACGTCCA**  
**GGAGCGCACCATCTTCTTCAAGGACGACGGCAACTACAAGACCCGCGCCGAGGTGAAGTT**  
**CGAGGGCGACACCCTGGTGAACCGCATCGAGCTGAAGGGCATCGACTTCAAGGAGGACGG**  
**CAACATCTTGGGGCACAAGCTGGAGTACAACATCAGCCACAACGTCTATATCACCGCC**  
**GACAAGCAGAAGAACGGCATCAAGGCCAACTTCAAGATCCGCCACAACATCGAGGACGG**  
**CAGCGTGCAGCTCGCCGACCACTACCAGCAGAACACCCCCATCGGCGACGGCCCCGTGCT**  
**GCTGCCCGACAACCACTACCTGAGCACCCAGTCCGCCCTGAGCAAAGACCCCAACGAGAA**  
**GCGCGATCACATGGTCCTGCTGGAGTTCGTGACCGCCGCGGGATC****GGATCC****GGCACCGGT**  
**GTAGGGCTACGTAAAGATGATGCTGAACTGACGGCTGCCTTCAATAAGGCGCTTGCGGAGC**  
**TGCGTCAGGACGGCACCTACGACAAGATGGCGAAAAAGTATTCGACTTTAATGTCTACGG**  
**TGACGGTGGCAGTGGAGGGAGCGGTGGAAGTGGCGGAAGCGCGCTACCGGAGACGGTAC**  
**GTATCGGAACCGATAACCACTACGCACCGTTCTCATCGAAAGATGCTAAAGGTGATTTTGT**  
**GGCTTTGATATCGATCTCGGTAACGAGATGTGCAAACGGATGCAGGTGAAATGTACCTGGG**  
**TTGCCAGTGACTTTGACGCGCTGATCCCCTCACTGAAAGCGAAAAAAATCGACGCTATTAT**  
**TTCGTCGCTTTCCATTACCGATAAACGTCAGCAGGAGATTGCCTTCTCCGACAAGCTGTACG**  
**CCGCAGATTCTCGTTTGATTGCGGCCAAAGGTTACCGATTACGCCAACGCTGGATTCACTG**  
**AAAGGTAAACATGTTGGTGTGCTGCAGGGATCAACCCAGGAAGCTTACGCTAACGAGACC**  
**TGGCGTAGTAAAGGCGTGGATGTGGTGGCCTATGCCAACCAGGATTTGGTCTATTCCGATCT**  
**GGCTGCAGGACGTCTGGATGCTGCGTTACAAGATGAAGTTGCTGCCAGCGAAGGATTCTC**  
**AAGCAACCTGCTGGTAAAGATTTGCCTTTGCT****GTCGAC****GAGCTGTT****CACCGGGGTGGTGCC**  
**CATCCTGGTCGAGCTGGACGGCGACGTAAACGGCCACAAGTTCAGCGTGTCCGGCGAGGG**  
**CGAGGGCGATGCCACCTACGGCAAGCTGACCCTGAAGTTCATCTGCACCACCGGCAAGCT**  
**GCCCGTGCCCTGGCCCCACCCTCGTGACCACCTTCGGCTACGGCCTGATGTGCTTCGCCCGCT**  
**ACCCCGACCACATGAAGCAGCACGACTTCTTCAAGTCCGCCATGCCCGAAGGCTACGTCCA**  
**GGAGCGCACCATCTTCTTCAAGGACGACGGCAACTACAAGACCCGCGCCGAGGTGAAGTT**  
**CGAGGGCGACACCCTGGTGAACCGCATCGAGCTGAAGGGCATCGACTTCAAGGAGGACGG**  
**CAACATCTTGGGGCACAAGCTGGAGTACAACATCAGCCACAACGTCTATATCATGGC**  
**CGACAAGCAGAAGAACGGCATCAAGGTGAACTTCAAGATCCGCCACAACATCGAGGACG**  
**GCAGCGTGCAGCTCGCCGACCACTACCAGCAGAACACCCCCATCGGCGACGGCCCCGTGC**  
**TGCTGCCCGACAACCACTACCTGAGCTACCAGTCCGCCCTGAGCAAAGACCCCAACGAGA**  
**AGCGCGATCACATGGTCCTGCTGGAGTTCGTGACCGCCGCGGGATCACTCTCGGCATGGA**  
**CGAGCTGTACAAGTAA**

## 6.2. Protein Sequence

MRGSHHHHHHGMADTRIGVTIYKSAGMVSKGEELFTGVVPILVELDGDVNGHKFSVSGEG  
EGDATYGKLTCLKFICTTGKLPVPWPTLVTTLTWGVQCFSRYPDHMKQHDFFKSAMPEGYVQERT  
IFFKDDGNYKTRAEVKFEGDTLVNRIELKGIDFKEDGNILGHKLEYNYISHNVYITADKQKNGIKA  
NFKIRHNIEDGSVQLADHYQNTPIGDGPVLLPDNHYLSTQSALS KDPNEKRDMVLLEFVTAA  
GIGSGTGVGLRKDDAELTAAFNKALGELRQDGTYDKMAKKYFDNFVYGDGGSGGSGGSGGSAL  
PETVRIGTDTTYAPFSSKDAKGDFVGFIDLGNECMCKRMQVKCTWVASDFDALIPSLKAKKIDAI  
SSLSITDKRQQEIAFSKLYAADSRLIAAKGSPIQPTLDSLKGKHVGVVLQGSTQEAYANETWRSKG  
VDVVAYANQDLVYSDLAAGRLDAALQDEVAASEGFLKQPAAGKDFAFVDELFTGVVPILVELD  
GDVNGHKFSVSGEGEGDATYGKLTCLKFICTTGKLPVPWPTLVTTFTGYGLMCFARYPDHMKQHD  
FFKSAMPEGYVQERTIFFKDDGNYKTRAEVKFEGDTLVNRIELKGIDFKEDGNILGHKLEYNYNS  
HNVIYIMADKQKNGIKVNFKIRHNIEDGSVQLADHYQNTPIGDGPVLLPDNHYLSYQSALS KDP  
NEKRDMVLLEFVTAAAGITLGMDELYK

## 7. Table S1

**Table S1.** Overview of the binding parameters of the toolbox sensors.

| Sensor              | 00      | 0F        | 0R        | F0     | FF      | FR      | R0   | RF     | RR     |
|---------------------|---------|-----------|-----------|--------|---------|---------|------|--------|--------|
| R <sub>0</sub>      | 1.70    | 1.95      | 2.30      | 1.52   | 1.49    | 1.92    | 2.20 | 1.74   | 2.28   |
| R <sub>sat</sub>    | 2.16    | 2.25      | 2.37      | 2.59   | 2.61    | 2.52    | 2.20 | 3.20   | 2.87   |
| ΔR                  | 0.46    | 0.30      | 0.07      | 1.07   | 1.12    | 0.60    | --   | 1.46   | 0.59   |
| K <sub>d</sub> (μM) | 107 ± 5 | 4.7 ± 0.8 | 2.5 ± 0.3 | 67 ± 3 | 3 ± 0.3 | 3 ± 0.2 | --   | 81 ± 2 | 27 ± 4 |

Includes the FRET-ratios in the non-bound (R<sub>0</sub>) state and under saturating conditions (R<sub>sat</sub>), the sensitivity (ΔR), and the affinity (K<sub>d</sub>) of the sensor variants for lysine.

## 8. Table S2

**Table S2.** Overview of the binding parameters of the sensor prototype without additional linkers (00) in fresh medium and in culture supernatant (refers to Figure 5 in the main paper).

| Measurement System  | R <sub>0</sub> | R <sub>sat</sub> | ΔR   | K <sub>d</sub> | pH  |
|---------------------|----------------|------------------|------|----------------|-----|
| Fresh medium        | 0.63 ± 0.01    | 0.84 ± 0.01      | 0.21 | 0.37 ± 0.04 mM | 7.0 |
| Culture supernatant | 0.93 ± 0.02    | 1.08 ± 0.00      | 0.15 | 0.63 ± 0.14 mM | 7.5 |
| MOPS buffer         | 1.70 ± 0.00    | 2.16 ± 0.00      | 0.46 | 107 ± 5 μM     | 7.3 |

The mean values of three independent measurements are shown. They were performed directly in the Biolector® cultivation system in the Flowerplates® at 1000 rpm and 30 °C. For comparison, the values measured in MOPS buffer are shown.

## 9. Setup and Data Analysis of Sensor Application in Microscale Cultivation Experiments of a Lysine Producer

In the following tables, the setup and analysis of the sensor application for L-lysine estimation in microtiter cultivation is explained in detail. See Table S3 for information about the plate layout. The fluorescence signals were measured constantly at λ<sub>Ex</sub> = 430 ± 5 nm, λ<sub>Em</sub> = 468 ± 5 nm (ECFP) and λ<sub>Ex</sub> = 430 ± 5 nm, λ<sub>Em</sub> = 532 ± 5 nm (Citrine), so the fluorescence intensity in each sample was recorded before each sampling and was subtracted from the measured fluorescence intensity in the presence of the sensor protein (Table S4). With the calculated FRET-ratios of the standards, the mean values and standard deviations of the in-plate calibration respective calibration lines were calculated (Table S5, Figure S5). Based thereon, the respective lysine concentrations in the samples were estimated (Table S6).

## 10. Table S3

Table S3. Layout of the cultivation plate (Flowerplate®).

| Well     | 01                                                                                          | 02 | 03 | 04 | 05                                    | 06                                    | 07                                     | 08                                      |
|----------|---------------------------------------------------------------------------------------------|----|----|----|---------------------------------------|---------------------------------------|----------------------------------------|-----------------------------------------|
|          | Samples                                                                                     |    |    |    | Calibration Standards                 |                                       |                                        |                                         |
| A = 0 h  | cultivation of L-lysine producing <i>C. glutamicum</i><br>DM1933 (4 wells per sample point) |    |    |    | 0 mM L-lysine<br>in-plate calibration | 1 mM L-lysine<br>in-plate calibration | 10 mM L-lysine<br>in-plate calibration | 100 mM L-lysine<br>in-plate calibration |
| B = 4 h  |                                                                                             |    |    |    |                                       |                                       |                                        |                                         |
| C = 8 h  |                                                                                             |    |    |    |                                       |                                       |                                        |                                         |
| D = 12 h |                                                                                             |    |    |    |                                       |                                       |                                        |                                         |
| E = 16 h |                                                                                             |    |    |    |                                       |                                       |                                        |                                         |
| F = 20 h |                                                                                             |    |    |    |                                       |                                       |                                        |                                         |

Note: In each row, four wells were filled with culture broth (1–4) and four wells were filled with L-lysine standards for recalibration of the sensors depending on the pH and on the medium composition. Accordingly, the standards in row A were prepared with fresh media, whereas the standards in rows B–F contained increasing amounts of the culture supernatant, which was prepared before with a wild-type strain of *C. glutamicum*. The calibration standards in row B were prepared with a mixture of 67% fresh medium and 33% culture supernatant and accordingly the standards in row C were prepared with a mixture of 33% fresh medium and 67% culture supernatant, whereas rows D, E, and F were prepared with 100% culture supernatant.

## 11. Table S4

Table S4. Raw data of the fluorescence measurement of the yellow and blue channel at sampling times. In-plate calibration data is shaded in gray, respectively.

| Well Number | ECFP-Signal ( $\lambda_{Ex} = 430 \pm 5 \text{ nm}$ , $\lambda_{Em} = 468 \pm 5 \text{ nm}$ ) |                                         |                                                                                                | Citrine-Signal ( $\lambda_{Ex} = 430 \pm 5 \text{ nm}$ , $\lambda_{Em} = 532 \pm 5 \text{ nm}$ ) |                                         |                                                                                                | FRET-Ratio =<br>Citrine/ECFP |
|-------------|-----------------------------------------------------------------------------------------------|-----------------------------------------|------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|-----------------------------------------|------------------------------------------------------------------------------------------------|------------------------------|
|             | Raw Data<br>Before<br>Biosensor<br>Addition                                                   | Raw Data After<br>Biosensor<br>Addition | Fluorescence Signal of the Sensor<br>Minus Background Fluorescence<br>of the Cultivation Broth | Raw Data<br>Before<br>Biosensor<br>Addition                                                      | Raw Data After<br>Biosensor<br>Addition | Fluorescence Signal of the Sensor<br>Minus Background Fluorescence<br>of the Cultivation Broth |                              |
| 1           | 0.41                                                                                          | 11.38                                   | 10.97                                                                                          | 0.20                                                                                             | 7.64                                    | 7.44                                                                                           | 0.68                         |
| A02         | 0.35                                                                                          | 11.11                                   | 10.75                                                                                          | 0.17                                                                                             | 7.61                                    | 7.44                                                                                           | 0.69                         |

|     |      |       |              |      |       |              |             |
|-----|------|-------|--------------|------|-------|--------------|-------------|
| A03 | 0.35 | 11.04 | <b>10.68</b> | 0.18 | 7.55  | <b>7.37</b>  | <b>0.69</b> |
| A04 | 0.38 | 11.11 | <b>10.73</b> | 0.18 | 7.52  | <b>7.34</b>  | <b>0.68</b> |
| A05 | 0.37 | 9.82  | <b>9.45</b>  | 0.17 | 6.57  | <b>6.40</b>  | <b>0.68</b> |
| A06 | 0.38 | 9.53  | <b>9.15</b>  | 0.17 | 6.17  | <b>6.00</b>  | <b>0.66</b> |
| A07 | 0.41 | 11.38 | <b>10.97</b> | 0.19 | 9.50  | <b>9.31</b>  | <b>0.85</b> |
| A08 | 0.65 | 13.92 | <b>13.27</b> | 0.29 | 12.08 | <b>11.79</b> | <b>0.89</b> |
| B01 | 0.55 | 12.39 | <b>11.85</b> | 0.40 | 8.38  | <b>7.98</b>  | <b>0.67</b> |
| B02 | 0.52 | 11.98 | <b>11.46</b> | 0.37 | 8.10  | <b>7.73</b>  | <b>0.67</b> |
| B03 | 0.55 | 12.15 | <b>11.60</b> | 0.38 | 8.34  | <b>7.96</b>  | <b>0.69</b> |
| B04 | 0.55 | 12.13 | <b>11.59</b> | 0.40 | 8.28  | <b>7.88</b>  | <b>0.68</b> |
| B05 | 0.66 | 14.13 | <b>13.47</b> | 0.77 | 9.88  | <b>9.11</b>  | <b>0.68</b> |
| B06 | 0.69 | 14.05 | <b>13.37</b> | 0.79 | 10.07 | <b>9.28</b>  | <b>0.69</b> |
| B07 | 0.72 | 14.97 | <b>14.26</b> | 0.81 | 12.14 | <b>11.33</b> | <b>0.79</b> |
| B08 | 0.99 | 16.84 | <b>15.85</b> | 0.89 | 14.58 | <b>13.69</b> | <b>0.86</b> |
| C01 | 0.95 | 12.09 | <b>11.14</b> | 0.99 | 9.08  | <b>8.09</b>  | <b>0.73</b> |
| C02 | 0.95 | 11.92 | <b>10.97</b> | 1.00 | 9.01  | <b>8.01</b>  | <b>0.73</b> |
| C03 | 0.98 | 12.32 | <b>11.34</b> | 1.00 | 9.26  | <b>8.26</b>  | <b>0.73</b> |
| C04 | 0.99 | 12.41 | <b>11.42</b> | 1.01 | 9.43  | <b>8.42</b>  | <b>0.74</b> |
| C05 | 0.95 | 16.81 | <b>15.86</b> | 1.20 | 11.66 | <b>10.46</b> | <b>0.66</b> |
| C06 | 0.98 | 16.66 | <b>15.69</b> | 1.21 | 12.01 | <b>10.80</b> | <b>0.69</b> |
| C07 | 1.06 | 17.08 | <b>16.02</b> | 1.24 | 13.81 | <b>12.57</b> | <b>0.78</b> |
| C08 | 1.35 | 19.05 | <b>17.70</b> | 1.33 | 16.81 | <b>15.48</b> | <b>0.87</b> |
| D01 | 1.46 | 12.50 | <b>11.04</b> | 1.79 | 10.26 | <b>8.47</b>  | <b>0.77</b> |
| D02 | 1.46 | 12.44 | <b>10.98</b> | 1.80 | 10.31 | <b>8.51</b>  | <b>0.78</b> |
| D03 | 1.49 | 12.52 | <b>11.02</b> | 1.83 | 10.49 | <b>8.66</b>  | <b>0.79</b> |

|     |      |       |       |      |       |       |      |
|-----|------|-------|-------|------|-------|-------|------|
| D04 | 1.52 | 12.72 | 11.20 | 1.84 | 10.65 | 8.81  | 0.79 |
| D05 | 1.38 | 20.36 | 18.98 | 1.77 | 14.45 | 12.68 | 0.67 |
| D06 | 1.36 | 20.14 | 18.78 | 1.76 | 14.65 | 12.89 | 0.69 |
| D07 | 1.39 | 20.20 | 18.81 | 1.77 | 16.09 | 14.32 | 0.76 |
| D08 | 1.62 | 21.03 | 19.41 | 1.87 | 18.27 | 16.40 | 0.84 |
| E01 | 1.36 | 13.17 | 11.80 | 1.96 | 11.00 | 9.04  | 0.77 |
| E02 | 1.35 | 13.07 | 11.72 | 1.95 | 11.05 | 9.10  | 0.78 |
| E03 | 1.40 | 13.22 | 11.82 | 1.97 | 11.20 | 9.23  | 0.78 |
| E04 | 1.40 | 13.32 | 11.92 | 2.00 | 11.31 | 9.31  | 0.78 |
| E05 | 1.42 | 20.78 | 19.36 | 1.82 | 14.38 | 12.56 | 0.65 |
| E06 | 1.40 | 20.67 | 19.26 | 1.78 | 14.81 | 13.03 | 0.68 |
| E07 | 1.43 | 20.12 | 18.69 | 1.80 | 16.33 | 14.53 | 0.78 |
| E08 | 1.69 | 21.41 | 19.72 | 1.89 | 18.43 | 16.54 | 0.84 |
| F01 | 1.43 | 13.28 | 11.85 | 2.07 | 11.01 | 8.94  | 0.75 |
| F02 | 1.42 | 13.37 | 11.95 | 2.06 | 11.15 | 9.09  | 0.76 |
| F03 | 1.42 | 13.38 | 11.96 | 2.10 | 11.25 | 9.15  | 0.77 |
| F04 | 1.48 | 13.34 | 11.86 | 2.14 | 11.26 | 9.12  | 0.77 |
| F05 | 1.42 | 20.43 | 19.01 | 1.79 | 14.34 | 12.55 | 0.66 |
| F06 | 1.43 | 20.67 | 19.24 | 1.81 | 14.63 | 12.82 | 0.67 |
| F07 | 1.39 | 19.80 | 18.41 | 1.76 | 15.52 | 13.76 | 0.75 |
| F08 | 1.66 | 20.70 | 19.04 | 1.87 | 18.01 | 16.14 | 0.85 |

## 12. Table S5

**Table S5.** FRET-ratios of the in-plate calibration with mean values and standard deviation.

| FRET-Ratios of the L-Lysine Standards |      |      |      |      |      |      | Mean Value  | Standard Deviation | Deviation in Percent |
|---------------------------------------|------|------|------|------|------|------|-------------|--------------------|----------------------|
| Sampling Time                         | 0 h  | 4 h  | 8 h  | 12 h | 16 h | 20 h |             |                    |                      |
| 0 mM L-lysine                         | 0.68 | 0.68 | 0.66 | 0.67 | 0.65 | 0.66 | <b>0.67</b> | 0.01               | 1.5                  |
| 1 mM L-lysine                         | 0.66 | 0.69 | 0.69 | 0.69 | 0.68 | 0.67 | <b>0.68</b> | 0.01               | 2.0                  |
| 10 mM L-lysine                        | 0.85 | 0.79 | 0.78 | 0.76 | 0.78 | 0.75 | <b>0.79</b> | 0.03               | 4.1                  |
| 100 mM L-lysine                       | 0.89 | 0.86 | 0.87 | 0.84 | 0.84 | 0.85 | <b>0.86</b> | 0.02               | 2.1                  |

Note: the saturation concentration of the sensor with L-lysine under the given conditions is  $> 10$  mM but  $< 100$  mM L-lysine and could not be resolved in the three-point calibration. Therefore, the calibration curves (Figure. S5) apparently increase continuously until 100 mM.

### 13. Figure S5




**Figure S5.** In-plate measured calibration curves.

**14. Table S6**

**Table S6.** Analysis of the culture broth using the sensor prototype.

| Sampling Time                               | 0 h  | 4 h  | 8 h  | 12 h | 16 h | 20 h |
|---------------------------------------------|------|------|------|------|------|------|
| FRET-ratios of the culture broth            |      |      |      |      |      |      |
| Well 1                                      | 0.68 | 0.67 | 0.73 | 0.77 | 0.77 | 0.75 |
| Well 2                                      | 0.69 | 0.67 | 0.73 | 0.78 | 0.78 | 0.76 |
| Well 3                                      | 0.69 | 0.69 | 0.73 | 0.79 | 0.78 | 0.77 |
| Well 4                                      | 0.68 | 0.68 | 0.74 | 0.79 | 0.78 | 0.77 |
| Mean Value                                  | 0.69 | 0.68 | 0.73 | 0.78 | 0.78 | 0.76 |
| Standard Deviation (SD)                     | 0.01 | 0.00 | 0.00 | 0.01 | 0.01 | 0.01 |
| Regression curve: $y = a \times \ln(x) + b$ |      |      |      |      |      |      |

|                           |                  |                  |                   |                    |                   |                 |
|---------------------------|------------------|------------------|-------------------|--------------------|-------------------|-----------------|
| <b>a</b>                  | 0.0506           | 0.0368           | 0.0405            | 0.0344             | 0.0352            | 0.0394          |
| <b>b</b>                  | 0.6812           | 0.6996           | 0.6895            | 0.6849             | 0.6831            | 0.6632          |
| L-lysine in culture broth |                  |                  |                   |                    |                   |                 |
|                           | 0–1 ± 0.12<br>mM | 0–1 ± 0.01<br>mM | 3 mM ±<br>0.01 mM | 15 mM ±<br>0.02 mM | 15 mM ±<br>0.3 mM | 12 ± 0.04<br>mM |

## References

1. Kremers, G.J.; Goedhart, J.; van Munster, E.B.; Gadella, T.W., Jr. Cyan and yellow super fluorescent proteins with improved brightness, protein folding, and FRET Forster radius. *Biochemistry* **2006**, *45*, 6570–6580.
2. Griesbeck, O.; Baird, G.S.; Campbell, R.E.; Zacharias, D.A.; Tsien, R.Y. Reducing the environmental sensitivity of yellow fluorescent protein. *J. Biol. Chem.* **2001**, *276*, 29188–29194.